

INDIVIDUAL COMPETENCY TASKING LIST

Nephrology/Dialysis, Walter Reed Army Medical Center, Washington, D.C. 20307-5001

Person's Name: _____ Rank/Grade: _____

Assigned Work Area: Nephrology Service Nephrology/Dialysis Lab

Specialty Area: Activated Clotting Time Quality Control Clot Sensitivity/Calibration

Indicate (by checking either "YES" or "NO" in the columns below) whether the person named above is required to demonstrate competency on the tasking list below. When the competency test has been successfully achieved, enter the date and the initials of the Authorized Supervisor entering the result.

Competency Tasking List	YES	NO	Date Completed	Auth Init.
Check that power is on				
Depress start button twice to start timer				
Insert an empty Hemochron test tube into test well				
Rotate tube twice counterclockwise to verify light illumination				
Wait 20 seconds				
Rotate slowly the entire instrument counterclockwise onto its left side until instrument reaches 90degrees from normal horizontal position				
Listen for audible beep/note flashing time on display screen				
Accept/reject results				
Record results in QC log book				

List approved by: _____ Date: _____